



## IMMUNOGEN SUBMISSION FORM

Please submit the form via Email: [Inquiries@JosmanLLC.com](mailto:Inquiries@JosmanLLC.com)

Sponsor's Name:	Date:
Immunogen ID:	Buffer:
Total amount submitted:	Concentration:
Number of tubes:	

Species of Animals: Rat <input type="checkbox"/> Mouse <input type="checkbox"/>
Dosage to be used for initial immunization per animal:
Dosage to be used per boost per animal:

Immunogen: Solution <input type="checkbox"/> Suspension <input type="checkbox"/> Lyophilized <input type="checkbox"/> Other <input type="checkbox"/> specify:
If any remaining Immunogen at completion: Discard <input type="checkbox"/> Return <input type="checkbox"/>

This immunogen does not contain live pathogens:  (Sponsor's initials)

Any material of human origin has been tested by a F.D.A. approved method and found to be non-reactive for HIV-I antibody and Hepatitis B surface antigen (HBs Ag):  (Sponsor's initials)

Please check which Storage Conditions are to be used at Josman, LLC: 4° C  -20° C  -80° C  Other:

### Recommended Guidelines:

Host Species:	Dosage/Animal/Injection:	Buffer:	Volume/dose
Mice	0.025mg – 0.05mg	PBS or other non-toxic	< 0.125 mL
Rats	0.05mg – 0.1mg	PBS or other non-toxic	≤ 0.250 mL

Please aliquot immunogen in doses to be given per animal per injection. Each aliquot should contain enough doses for all animals on the project. For example: If there are 2 rabbits on the project, each aliquot should contain both doses.

Special instructions/comments:
--------------------------------

Telephone: (707)226-2115  
Email: [Inquiries@JosmanLLC.com](mailto:Inquiries@JosmanLLC.com)  
[www.JosmanLLC.com](http://www.JosmanLLC.com)